ANKLE BRACE PRO COMFORT

GB IE

ANKLE BRACE

Instructions for use

FI

NILKKATUKISIDE

Käyttöohje

SE

FOTLEDSBANDAGE

Bruksanvisning

DK

FODLEDSBANDAGER

Brugsanvisning

IAN 93485

GB IE FI SE DK

Congratulations on the purchase of your new product. You have chosen a high-quality article. The instructions for use are part of this product. They contain important information on safety and use. Please familiarise yourself with the instructions and safety information before using this product. Use the product only in perfect condition, as described and for the specified area of use. Please keep all packaging including this leaflet for as long as you use the brace and be sure to include them when handing on the product to another person.

For personal use only. Consult a doctor if:

- you are unsure how to use the product
- you have acute injuries
- unforeseen problems arise

Consult a doctor before wearing the brace for hours at a time and do not wear in combination with other medical products, e.g. creams or salves. Ensure good circulation when wearing the brace!

	Consult instructions for use
	Manufacturing date
	Batch designation
	Manufacturer

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
GB IE

Applying the brace

To reduce the tendency towards supination (an outward rolling of the foot) and pronation (an inward rolling of the foot)

Image 1: Place your heel on the heel ring

Image 2: First close the bottom and then the top Velcro fastener

With a tendency towards supination

Image 3: First take the outer stabilisation strap and pull it from the outside inwards across the arch of the foot (instep) and close the Velcro fastener, stretching it slightly.

Image 4: Then pull the inner stabilisation strap across the arch of the foot towards the ankle and close the Velcro fastener. The tighter the stabilisation strap is pulled, the more supination is reduced.

With a tendency towards pronation

Image 5: First take the inner stabilisation strap and pull it from the inside outwards across the arch of the foot (instep) and close the Velcro fastener, stretching it slightly.

Image 6: Then pull the outer stabilisation strap across the arch of the foot towards the inside of the ankle and close the Velcro fastener. The tighter the stabilisation strap is pulled, the more pronation is reduced.

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Horizonte Textil GmbH, Am Sandtorkai 23, D-20457 Hamburg