GB IE

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Congratulations! With your purchase you have decided on a high-quality product. Get to know the product before you start to use it. Carefully read the following instruction manual. Use the product only as described and only for the given areas of application. Keep these instructions safe. When passing the product on to a third party, always make sure that the documentation is included.

## Scope of delivery (figure A)

- 1 x Mobility trainer with resistance screw and display (including LR44 battery)
- 1 x Instruction manual

## Technical data

Dimensions: approx. 48 x 37 x 30cm

Weight: approx. 2.5kg

Material: plastic, steel

Maximum load capacity: 100kg

Equipment class: HC (for use in private households, precision C) / not suitable for therapeutic training

Speed-independent training device

## Correct use

The article is not designed for commercial use. This article enables you to train your arms and legs; fitted with a resistance screw in order to increase and maintain mobility. An integrated display shows training advances. Folding legs allow you to store the mobility trainer and assemble it again.

## Safety notes

**Important: read through these operating instructions carefully and always retain them for future reference.**

- The article is not suitable for therapeutic purposes.
- During the exercises sit upright and relaxed with a straight back.
- Exchange defective parts of the article immediately and do not use the article until the repair is complete.


- Please note that incorrect or excessive training may be detrimental to health.
- Place the article on a horizontal, level surface.
- Use an anti-slip underlay on smooth surfaces.
- Observe wear on parts. If the article no longer offers any step resistance then the drive belt is worn and the article needs to be replaced.
- There must be a free area of approx. 0,6m around the training area (figure H, I).

## **⚠ Danger to life!**

- Never allow children to play with the packaging material unsupervised. There is a risk of suffocation.
- The product is not a piece of climbing or play equipment. Make sure that persons, particularly children, do not stand on the product or pull themselves up on it. The product can tip over.
- The article is not designed to support the entire body weight.
- Children should only use the article under supervision.
- Consult your doctor before beginning use or if you experience any discomfort.

## **⚠ Risk of injury!**

- Do not place your full body weight on the article. It may only be used as described in the instructions.

- Ensure that the article is sufficiently stable before use.
- Place the article on a level surface.
- Be careful of your fingers when unfolding and folding. Fingers may become trapped.
- Ensure that the article is fully extended before use.
- Examine the article for any damage or wear before each use. Replace damages or worn parts on the article immediately as the function and safety could otherwise be impeded.

- The article may only be used for the intended purpose.
- Do not touch the metal parts on the article after use, particularly not the pedal arms as they heat up during use.

### **Warning notes - batteries!**

- Only use batteries from the same manufacturer and of the same type.
- Remove the battery when it is spent or if the product has not been used for a long period.
- Warning! Batteries may not be charged or re-activated by other means, may not be dismantled, thrown into the fire, or short-circuited.
- Always keep batteries out of the reach of children.
- Non re-chargeable batteries may not be charged.
- Re-chargeable batteries may only be charged under adult supervision.
- Re-chargeable batteries must be removed from the article before being charged.
- The connecting terminals may not be short-circuited.
- Treat a damaged or leaking battery with great care and dispose of it immediately in the correct way.
- If you come into contact with battery acid, then wash the affected area with water and soap. If battery acid gets in your eyes then rinse them with water and consult a doctor immediately.

## **Assembly**

### **Assembling the article**

1. Remove the packaging material.
2. Unfold both legs until you hear an audible click from the safety pins (figure B).
3. Put the article up in front of you and turn the resistance screw clockwise to increase resistance and counter-clockwise to decrease resistance (figure C).

**Note:** '+' and '-' are marked on the resistance screw to adjust the resistance. '+' stands for increase and '-' stands for reduction of the resistance.

### **Folding the article (figure D)**

1. Turn the article onto its side.
2. Pull on the key ring on the legs to release the safety pin and fold the legs inwards, one at a time.

### **Setting the pedal strap (figure E)**

The pedal straps are pre-mounted. You have the choice of setting the straps to one of two lengths.

1. Pull the outer strap end over the second pedal pins.
2. Choose the holes on the strap ends for the appropriate length and place them over the pedal pins.

### **Changing the battery**

**CAUTION! Follow these instructions in order to avoid mechanical and electrical damage.**

#### **Danger!**

One LR44 battery is included in the article. If the battery is no longer working then you can change it.

1. Lift the display on the side and remove it carefully (figure F). The battery is on the back of the display.
2. Carefully slide the spent battery out of the battery holder and exchange it for a new battery (figure G).

**Note:** ensure that the battery is completely inside the mounting.

3. Insert the cable connections and the display back into the casing.

**Note:** please ensure that the display is replaced correctly into the casing.

### **Training notes**

- Consult a doctor before using the article.
- Stop training as soon as you feel discomfort during training.
- Start slowly and increase training intensity as necessary.
- We recommend training the legs and arms every day for 30 minutes each.

- Wear comfortable clothing to use the article. However please ensure that trouser legs or sleeves are not too wide as they could be caught in the pedals.
- Do some easy stretches before each use, and start with slow, gentle steps. Increase slowly until you have reached the desired speed.
- Only increase resistance slowly.
- Slow the speed to end use.

## Use

Put the article in front of you ready for training. Start with easy training and resistance levels. Release the resistance screw. The tighter the resistance screw is turned, the more difficult it is to turn the pedals. Warm up before training by carrying out stretching and warm-up exercises. You should relax after each training session.

## Leg training (figure H)

1. Put the article on a level, stable and non-slip floor and sit in front of the article on a chair, relaxed and upright with a straight back.

**Note:** ensure that the chair is stable and comfortable.

2. Place your feet on the pedals in such a way that the pedal strap is over your foot.
3. 'Pedal' your feet as if on a bicycle.

## Arm training (figure I)

1. Put the article in front of you on a level, stable and non-slip surface, e.g. at table height, and sit on a chair, relaxed and upright with a straight back.
2. Place your palms on the pedals in such a way that the pedal strap is over your hand.
3. 'Pedal' your arms as if on a bicycle.

## Display function

The display has various modes that you can select by pressing the button:

- SCAN: This mode displays the various display modes alternately.
- TIME: This mode displays the duration of your current training session.
- CNT: This mode displays the rounds completed in the current training session.

- CAL: This mode displays the calories used during your current training session.

**Note:** calorie consumption is only an approximate reference value.

- R.P.M.: This mode displays the speed in rounds per minute.

**Note:** the display switches off automatically after approx. four minutes when not in use.

## Set value to zero

Press on the button for approximately four seconds to reset the value to 0.

## Cleaning and care

Clean the article with a damp cloth and a mild cleaning agent.

Caution! Do not allow moisture to penetrate the display. Fold the legs away as described in the section 'folding the article'. Store the article in a dry place.

Parts that are susceptible to wear: if the article is no longer offering any step resistance then the drive belt is worn. The article must be replaced.

## Information on disposal

Please dispose of packaging and product in an environmentally friendly way and according to the material type!

Dispose of the product through an authorised waste disposal company or your municipal authorities. Please observe the currently applicable regulations.

## Disposing of batteries

Batteries must not be disposed of along with normal domestic waste!

Return batteries to recycling facilities or the dealer who sold them.

Batteries are marked with a dustbin with a red cross through it and chemical pollutant symbols „Cd“ for cadmium, „Hg“ for mercury and „Pb“ for lead. Contact your local authority for further information if in doubt.

## 3 Years Warranty

The product was produced with great care and under constant supervision. You receive a three-year warranty for this product from the date of purchase. Please retain your receipt.

The warranty applies only to material and workmanship and does not apply to misuse or improper handling. Your statutory rights, especially the warranty rights, are not affected by this warranty.

With regard to complaints, please contact the following service hotline or contact us by e-mail.

Our service employees will advise as to the subsequent procedure as quickly as possible. We will be personally available to discuss the situation with you.

Any repairs under the warranty, statutory guarantees or through goodwill do not extend the warranty period. This also applies to replaced and repaired parts. Repairs after the warranty are subject to a charge

IAN: 96826

### **GB** Service Great Britain

Tel.: 0871 5000 720  
(£ 0.10/Min.)

E-Mail: [deltasport@lidl.co.uk](mailto:deltasport@lidl.co.uk)

### **IE** Service Ireland

Tel.: 1890 930 034  
(0,08 EUR/Min., (peak))  
(0,06 EUR/Min., (off peak))

E-Mail: [deltasport@lidl.ie](mailto:deltasport@lidl.ie)

You can also find spare parts for your product at:  
[www.delta-sport.com](http://www.delta-sport.com), category Service - Lidl Spare Parts  
Service

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