RÜCKEN-NACKEN-HEIZKISSEN / BACK & NECK HEAT PAD / COUSSIN CHAUFFANT DOS ET NUQUE SRNH 100 H5

DE AT CH

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FR BE

COUSSIN CHAUFFANT DOS ET NUQUE

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PL

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IMPORTANT INSTRUCTIONS – KEEP FOR FUTURE REFERENCE!

Back & neck heat pad

● Introduction

This quick reference guide is a fixed part of the operating instructions. Keep both this and the operating instructions in a safe place. When passing this product on to third parties, please be sure to include all documentation. Read the operating instructions before use, and pay particular attention to the safety instructions included in the operating instructions.

Scan the QR code, or download the complete operating instructions at www.lidl-service.com.

● Intended use

This heating pad is intended for applying heat to the human body. It can be used to apply heat to the back & neck. This heating pad is not suitable for use in hospitals or for any other commercial use. In particular it must not be used for babies or infants, persons who are insensitive to heat, helpless persons,

or animals. Any use other than that described above or modification of the heating pad are not permissible and may lead to injury and/or damage to the heating pad. The manufacturer assumes no liability for damage arising out of improper use.

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Important safety instructions

WARNING! Failure to observe the following instructions can cause personal injury or property damage (electric shock, burns, fire). The following safety instructions and hazard warnings are not only intended to protect your own health and the health of others, but also to prevent damage to the heating pad. Please observe these safety instructions and pass on this operating manual together with the heating pad.

SCALDING HAZARD!

This heating pad may not be used by persons insensitive to heat and other persons in need

- of protection who do not respond to overheating (e.g. diabetics, persons with skin changes caused by illness or scarred skin areas in the application area, after using analgesics or alcohol). There is a risk of receiving skin burns.
-  **WARNING!** This heating pad may not be used for very young children (0–3 years), as they are unable to react to overheating. There is a risk of receiving skin burns.
 - The heating pad may be used by children over 3 years and under 8 years on the lowest setting when supervised. Burns may otherwise occur.
 - This heating pad may be used by older children age 8 and up, and persons with diminished physical, sensory or mental capacities or lacking experience and knowledge when supervised and instructed on the safe use of the heating pad and understand the associated risks.
 - Never allow children to play with the heating pad. There is a risk of injury.
 - Children must not perform clean or maintain the product without supervision. There is a risk of injury.
 - This heating pad is not intended for hospital or commercial use. This poses a risk of injuries and/or damage to the heating pad.
 -  **WARNING!** Do not insert needles or pointed objects into the heating pad. There is a risk of electric shock.
 -  **WARNING!** Do not switch on the heating pad whilst folded or pushed together. The heating pad may otherwise be damaged.
 - Do not use when wet! Otherwise there is a risk of electric shock.
 - This heating pad may only be operated together with the control unit specified on the label. The heating pad may otherwise be damaged.
 - The electric and magnetic fields emitted by this heating pad could potentially interfere with the functionality of your pacemaker. However, the levels are far below the permissible limits: electric field strength: max. 5000V/m, magnetic field strength: max. 80 A/m, magnetic flux density: max. 0.1 milli-Tesla. Therefore please consult your physician and the manufacturer of your pacemaker before using this heating pad.

- Do not pull, twist or kink the cables. The heating pad may otherwise be damaged.
- The cable and the switch of the heating pad may pose risk of catching, strangulation, tripping or being stepped on if positioned incorrectly. The user must ensure excess cable and cables in general are positioned safely.
- Frequently check the heating pad for wear or damage. If signs thereof are detected, the heating pad was used improperly or will no longer heat up, it must first be inspected by the manufacturer before being switched on again.
- If the power cord of the heating pad is damaged, it must be replaced by the manufacturer or its service department, or a similarly qualified person, to avoid hazards.
- With the heating pad switched on, never
 - place objects on it (e.g. suitcases or clothes baskets),
 - place a heating source such as a hot-water bottle, heating pad or similar items on it. The heating pad may otherwise be damaged.
- The electronic components inside the control unit will heat up whilst

operating the heating pad. For this reason, the control unit must not be covered or be lying on the heating pad during operation. This poses a risk of injuries and / or damage to the heating pad.

- Always observe the information on start-up, cleaning and care, and storage.
- If you have additional questions on using our products, please contact our Customer Service Department.

● Start-up

Note: The heating pad may smell of plastic the first time you use it, but this will disappear after a short period of time.

- To begin, first connect the control unit **1** to the heating pad with the plug-in connector **3** (see fig. B).
- Position the heating pad on your back so that the fastener for the neck section is in the throat area.
- Then close the fastener.
- Then adjust the length of the belt around your stomach and fit the two ends of the connector together (see fig. C) to fasten the belt.
- To open the belt, press the two catches on the connector at the same time (see fig. C).

● Switching on / off / Selecting the temperature

Switching on:

- Set the control unit **1** to 1, 2, 3, 4, 5 or 6 to switch on the heating pad.

Note: The function display **2** illuminates after

the heating pad is switched on.

Note: This heating pad features rapid heating, which allows it to heat up within the first 10 minutes of operation.

Switching off:

- Set the control unit **1** to 0 to switch off the heating pad.

Note: The operating display **2** will switch off.

Selecting the temperature:

- Just before using the heating pad, set it to the highest temperature setting (setting 6). This will cause the heating pad to warm up quickly.
- When using the heating pad for several hours, set it to the lowest temperature setting (setting 1).

⚠ WARNING! When operating the heating pad for several hours, we commend setting the control unit **1** to the lowest setting to prevent overheating and potential burns to the part of the body heat is applied to.

● Storage

-  **RISK OF DAMAGE TO THE PRODUCT OR PROPERTY!** Let the heating pad cool before putting it away. The heating pad may otherwise be damaged.
-  **RISK OF DAMAGE TO THE PRODUCT OR PROPERTY!** To prevent the heating pad from becoming sharply creased, do not store other objects on top of the heating pad. The heating pad may otherwise be damaged.

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NU Service GmbH